AQUA LAB

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13. Literatura

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- Bousquet-Ricard, M., G. Qualyle, T. Pharm, and J.C. Cheftel. (1980). Comparative study of three methods of determining water activity in intermediate moisture foods. *Lebensmittel- Wissenschaft und-Technologie*. 13:169-173.
- Cazier, J.B., and V. Gekas. 2001. Water activity and its prediction: a review. *International Journal of Food properties* 4(1):35-43.
- Chirife, J., G. Favetto, C. Ferro-Fontan, and S. Resnik. (1983). The water activity of standard saturated salt solutions in the range of intermediate moisture foods. *Lebensmittel-Wissenschaft und-Technologie*. 16:36-38.
- Duckworth, R. (1975). *Water Relations of Foods*. Academic Press, New York.
- Gomez-Diaz, R. (1992). Water activity in foods: Determination methods. *Alimentaria*. 29:77-82. Gómez, R. and Fernandez-Salguero J. (1992). Water activity and chemical composition of some food emulsions. *Food Chemistry*. 45:91-93.
- Gómez, R. and Fernandez-Salguero J. (1992). Water activity and chemical composition of some food emulsions. *Food Chemistry*. 45:91-93.
- Greenspan, L. (1977). Humidity fixed points of binary saturated aqueous solutions. *Journal of Research of the National Bureau of Standards A.Physics and Chemistry*. 81A:89-96.
- Karmas, E. (1981). Measurement of moisture content. *Cereal Foods World*. 26(7):332-334.
- Kitic, D., D.C. Pereira-Jardim, G. Favetto, S.L. Resnik, and J. Chirife. (1986). Theoretical prediction of the water activity of standard saturated salt solutions at various temperatures. *Journal of Food Science*. 51:1037-1042.
- Labuza, T.P. and R. Contreras-Medellin. (1981). Prediction of moisture protection requirements for foods. *Cereal Foods World*. 26(7):335-343.
- Labuza, T.P., K. Acott, S.R. Tatini, and R.Y. Lee. (1976). Water activity determination: A collaborative study of different methods. *Journal of Food Science*. 41:910-917.
- Marcolli, C., and Th . Peter. 2005. Water activity in polyol/water systems: new UNIFAC parameterization. *Atmospheric Chemistry and Physics* 5:1545-1555.
- Ninni, L., M.S. Camargo, and A.J.A. Meirelles. 2000. Water activity in polyol systems. *Journal of Chemical and Engineering Data* 45:654-660.
- Prior, B.A. (1979). Measurement of water activity in foods: A review. *Journal of Food Protection*. 42(8):668-674.
- Rahman, M.S. and S.S. Sablani. 2001. Measurement of water activity by electronic sensors. P. A2.5.1-A2.5.4 In R.E.Wrolstad (ed.) *Current Protocols In Food Analytical Chemistry*. John Wiley & Sons, Inc., New York.
- Rahman, M.S., S.S. Sablani, N. Guizani, T.P. Labuza, and P.P. Lewicki. 2001. Direct manometric determination of vapor pressure. P. A2.4.1-A2.4.6. In R.E. Wrolstad (ed.) *Current Protocols In Food Analytical Chemistry*. John Wiley & Sons, Inc., New York.

-
- Reid, D.S., A.J. Fontana, M.S. Rahman, S.S. Sablani, T.P. Labuza, N. Guizani, and P.P. Lewicki. 2001. Vapor pressure measurements of water p. A2.1.1-A2.5.4. In R.E. Wrolstad (ed.) *Current Protocols In Food Analytical Chemistry*. John Wiley & Sons, Inc., New York.
- Reid, D.S. (1976). Water activity concepts in intermediate moisture foods. In: *Intermediate Moisture Foods*. Davies, R., G.G. Birch, and K.J. Parker (ed.) Applied Science Publishers, London. pp. 54-65.
- Richard, J. and T.P. Labuza. (1990). Rapid determination of the water activity of some reference solutions, culture media and cheese using a new dew point apparatus. *Sciences des Aliments*. 10:57-64.
- Roa, V. and M.S. Tapia de Daza. (1991). Evaluation of water activity measurements with a dew point electronic humidity meter. *Lebensmittel- Wissenschaft und- Technologie*. 24(3):208-213.
- Rodel, W. 2001. Water activity and its measurement in food. P. 453-483. In E. Kress-Rogers, and C.B. Brimelow (ed.) *Instrumentation and sensors for the food industry*. CRC Press LLC, Boca Raton, FL.
- Roos, K.D. (1975). Estimation of water activity in intermediate moisture foods. *Food Technology*. 29:26-30.
- Scott, V.N. and D.T. Bernard. (1983). Influence of temperature on the measurement of water activity of food and salt systems. *Journal of Food Science*. 48:552-554.
- Snively, M.J., J.C. Price, and H.W. Jun. (1990). A comparison of three equilibrium relative humidity measuring devices. *Drug Development and Industrial Pharmacy*. 16(8):1399-1409.
- Stamp, J.A., S. Linscott, C. Lomauro, and T.P. Labuza. (1984). Measurement of water activity of salt solutions and foods by several electronic methods as compared to direct vapor pressure measurement. *Journal of Food Science*. 49:1139- 1142.
- Stoloff, L. (1978). Calibration of water activity measuring instruments and devices: Collaborative study. *Journal of Association of Official Analytical Chemists*. 61:1166-1178.
- Troller, J.A. (1983). Methods to measure water activity. *Journal of Food Protection*. 46:129-134.
- Troller, J.A. and J.H.B. Christian. (1978). *Water Activity and Food*. Academic Press, New York.
- Troller, J.A. and V.N. Scott. (1992). Measurement of water activity (A_w) and acidity. In: *Compendium of Methods for the Microbiological Examination of Foods*. Vanderzant, C. and D.F. Splittstoesser (ed.) American Public Health Association, Washington, D.C. pp. 135-151.
- van den Berg, C. (1985). Water activity. In: *Concentration and Drying of Foods*. MacCarthy, D. (ed.) Elsevier, London. pp. 11-35.
- van den Berg, C. (1991). Food-water relations: Progress and integration, comments and thoughts. In: *Water Relationships in Foods*. Levine, H. and L. Slade (ed.) Plenum Press, New York.
- van den Berg, C. and Bruin. (1981). Water activity and its estimation in food systems: Theoretical aspects. In: *Water Activity: Influences on Food Quality*. Rockland, L.B. and G.F. Stewart (ed.) Academic Press, New York. pp. 1-61.
-

- Vega, M.H. and G.V. Barbosa-Canovas. (1994). Prediction of water activity in food systems: A review on theoretical models. *Revista Espanola De Ciencia Y Tecnologia De Alimentos*. 34:368- 388.
- Vega, M.H., B. Romanach, and G.V. Barbosa-Canovas. (1994). Prediction of water activity in food systems: A computer program for predicting water activity in multicomponent foods. *Revista Espanola De Ciencia Y Tecnologia De Alimentos*. 34:427-440.
- Vos, P.T. and T.P. Labuza. (1974). Technique for measurements of water activity in the high a w range. *Journal of Agricultural and Food Chemistry*. 22:326-327.
- Voysey, P. (1993). An evaluation of the AquaLab CX-2 system for measuring water activity. *Digest, Microbiology Section*. 124

13.2. Kvalita a bezpečnost potravin

- Bei, Z.H., and R.-M.J. Nout. 2000. Effects of temperature, water activity and gas atmosphere on mycelial growth of tempe fungi *Rhizopus microsporus* var. *microsporus* and *R. microsporus* var. *oligosporus*. *World Journal of Microbiology and Biotechnology* 16:853-858.
- Beuchat, L.R. 1981. Microbial stability as affected by water activity. *Cereal Foods World* 26:345-349.
- Brandt, L. (1996). Bound for success. Controlling water activity gives technologists the edge in developing safe, shelf-stable foods. *Food Formulating*. September:41-48.
- Chirife, J. and B.M. Del-Pilar. (1994). Water activity, glass transition and microbial stability in concentrated/semimoist food systems. *Journal of Food Science*. 59:921-927.
- Chirife, J. and M.P. Buera. (1995). A critical review of some non-equilibrium situations and glass transitions on water activity values of foods in the microbiological growth range. *Journal of Food Engineering*. 25:531-552.
- Chirife, J. and M.P. Buera. (1996). Water activity, water glass dynamics, and the control of microbiological growth in foods. *Critical Reviews in Food Science and Nutrition*. 36(5):465-513.
- Farber, J.M., F. Coates, and E. Daley. 1992. Minimum water activity requirements for the growth of *Listeria monocytogenes*. *Lett Appl Microbiol* 15:103-105.
- Franks, F. (1982). Water activity as a measure of biological viability and quality control. *Cereal Foods World*. 27(9):403-407.
- Franks, F. (1991). Water activity: a credible measure of food safety and quality? *Trends in Food Science and Technology*. March:68-72.
- Garcia de Fernando, G.D., O. Diaz, M. Fernandez, and J.A. Ordonez. 1992. Changes in water activity of selected solid culture media throughout incubation. *Food Microbiology* 9:77-82.
- Gibson, A.M., J. Baranyi, J.I. Pitt, M.J. Eyles, and T.A. Roberts. 1994. Predicting fungal growth: The effect of water activity on *Aspergillus flavus* and related species. *International Journal of Food Microbiology* 23:419-431.
- Goaleni, N., J.E. Smith, J. Lacey, and G. Gettinby. 1997. Effects of temperature, water activity, and incubation time on production of aflatoxins and cyclopiazonic acid by

- an isolate of *Aspergillus flavus* in surface agar culture. *Appl Environ Microbiol* 63:1048-1053.
- Hardman, T.M. (1988). *Water and Food Quality*. Elsevier Press, London.
- Hocking, A.D., and B.F. Miscamble. 1995. Water relations of some Zygomycetes isolated from food. *Mycological Research* 99:1113-1118.
- Hocking, A.D., B.F. Miscamble, and J.I. Pitt. 1994. Water relations of *Alternaria alternata*, *Cladosporium cladosporioides*, *Cladosporium sphaerospermum*, *Curvularia lunata* and *Curvularia pallescens*. *Mycological Research* 98:91-94.
- Houtsma, P.C., A. Heuvelink, J. Dufrenne, and S. Notermans. 1994. Effect of sodium lactate on toxin production, spore germination and heat resistance of proteolytic *Clostridium botulinum* strains. *Journal of Food Protection* 57:327-330.
- Kress-Rogers, E. (1993). Food quality measurement. *Food Industry News*. September:23-26.
- Kuntz, L.A. 1992. Keeping microorganisms in control. *Food Product Design* August:44-51.
- Levine, H. and L. Slade. (1991). *Water Relationships in Foods*. Plenum Press, New York.
- Li, K.Y., and J.A. Torres. 1993. Water activity relationships for selected mesophiles and psychrotrophs at refrigeration temperature. *Journal of Food Protection* 56:612-615.
- Li, K.Y., and J.A. Torres. 1993. Water activity relationships for selected mesophiles and psychrotrophs at refrigeration temperature, *Journal of Food Protection* 56:612-615.
- Lopez-Malo, A., S. Guerrero, and S.M. Alzamora. 2000. Probabilistic modeling of *Saccharomyces cerevisiae* inhibition under the effects of water activity, pH, and potassium
- Mannheim, C.H., J.X. Liu, and S.G. Gilbert. (1994). Control of water in foods during storage. *Journal of Food Engineering*. 22:509-532.
- Marauska, M., A. Vigants, A. Klinec, D. Upite, E. Kaminska, and M. Bekers. 1996. Influence of water activity and medium osmolality on the growth and acid production of *Lactobacillus casei* var. *lactosus*. *Proceedings of the Latvian Academy of Sciences Section B Natural Exact and Applied Sciences* 50:144-146.
- Masana, M.O., and J. Baranyi. 2000. Growth/no growth interface of *Brochothrix thermosphacta* as a function of pH and water activity. *Food Microbiology* 17:485-858.
- Mattick, K. L., F. Jorgensen, J.D. Legan, M.B. Cole, J. Porter, H.M. Lappin-Scott, and T.J. Humphrey. 2000. Survival and filamentation of *Salmonella enterica* serovar Enteritidis PT4 and *Salmonella enterica* serovar Typhimurium DT104 at low water activity. *Appl Environ Microbiol* 66:1274-1279.
- Mattick, K.L., F. Jorgensen, J.D. Legan, H.M. Lappin-Scott, and T.J. Humphrey. 2000. Habituation of *Salmonella* spp. at reduced water activity and its effect on heat tolerance. *Appl Environ Microbiol* 66:4921-4925.
- Mattick, K.L., F. Jorgensen, J.D. Legan, H.M. Lappin-Scott, and T.J. Humphrey. 2001. Improving recovery of *Salmonella enterica* Sero var Typhimurium DT104 cells injured by heating at different water activity values. *Journal of Food Protection* 64:1472-1476.

- McMeekin, T.A. and T. Ross. (1996). Shelf life prediction: Status and future possibilities. *International Journal of Food Microbiology* 33:65-83.
- Miller, A.J. 1992. Combined water activity and solute effects on growth and survival of *Listeria monocytogenes*. *Journal of Food Protection* 55:414-418.
- Nakajo, M., and Y. Moriyama. 1993. Effect of pH and water activity on heat resistance of spores of *Bacillus coagulans*. *Journal of the Japanese Society for Food Science and Technology* 40:268-271.
- Nelson, K.A. and T.P. Labuza. (1994). Water activity and food polymer science: Implications of state on arrhenius and WLF models in predicting shelf life. *Journal of Food Engineering* 22:271-289.
- Nesci, A., M. Rodrigues, and M. Etcheverry. 2003. Control of *Aspergillus* growth and aflatoxin production using antioxidants at different conditions of water activity and pH. *Journal of Applied Microbiology* 95:279-287.
- Nolan, D.A., D.C. Chamblin, and J.A. Troller. 1992. Minimal water activity levels for growth and survival of *Listeria monocytogenes* and *Listeria innocua*. *International Journal of Food Microbiology* 16:323-335.
- Noorlidah, A., A. Nawawi, and I. Othman. 2000. Fungal spoilage of starch-based foods in relation to its water activity (aw). *Journal of Stored Products Research* 36:47-54.
- Park, C.M., and L.R. Beuchat. 2000. Survival of *Escherichia coli* O157:H7 in potato starch as affected by water activity, pH and temperature. *Lett Appl Microbiol* 31(5):364-367.
- Petersson, S., and J. Schnuerer. 1995. Biocontrol of mold growth in high-moisture wheat stored under airtight conditions by *Pichia anomala*, *Pichia guilliermondii*, and *Saccharomyces cerevisiae*. *Appl Environ Microbiol* 61:1027-1032.
- Pitt, J.I., and B.F. Mischamble. 1995. Water relations of *Aspergillus flavus* and closely related species. *Journal of Food Protection* 58:86-90.
- Plaza, P., J. Usall, N. Teixido, and I. Vinas. 2003. Effect of water activity and temperature on germination and growth of *Penicillium digitatum*, *P. italicum* and *Geotrichum candidum*. *Journal of Applied Microbiology* 94:549-554.
- Plaza, P., J. Usall, N. Teixido, and I. Vinas. 2003. Effect of water activity and temperature on germination and growth of *Penicillium digitatum*, *P. italicum* and *Geotrichum candidum*. *Journal of Applied Microbiology* 94:549-554.
- Quintavalla, S., and G. Parolari. 1993. Effects of temperature, water activity and pH on the growth of *Bacillus* cells and spore: A response surface methodology study. *International Journal of Food Microbiology* 19:207-216.
- Rockland, L.B., and G.F. Stewart. 1981. Water activity: Influences on food quality. Academic Press, New York.
- Rockland, L.B. and S.K. Nishi. (1980). Influence of water activity on food product quality and stability. *Food Technology* 34:42-59.
- Saad, R.R. 1992. Effect of water activity on growth and lipids of xerophilic fungi, *Aspergillus repens* and *Aspergillus amstelodami*. *Zentralblatt Fuer Mikrobiologie* 147:61-64.
- Salter, M.A., D.A. Ratkowsky, T. Ross, and T.A. McMeekin. 2000. Modelling the combined temperature and salt (NaCl) limits for growth of a pathogenic *Escherichia*

coli strain using nonlinear logistic regression. *International Journal of Food Microbiology* 61:159-167.

- Santos, J., T.M. Lopez-Diaz, M.C. Garcia-Lopez, M.C. Garcia-Fernandez, and A. Otero. 1994. Minimum water activity for the growth of *Aeromonas hydrophila* as affected by strain, temperature and humectant. *Lett Appl Microbiol* 19:76-78.
- Sautour, M., A. Rouget, P. Dantigny, C. Divies, and M. Bennsoussan. 2001. Prediction of conidial germination of *Penicillium chrysogenum* as influenced by temperature, water activity and pH. *Lett Appl Microbiol* 32:131-134.
- Seow, C.C., T.T. Teng, and C.H. Quah. (1988). *Food Preservation by Moisture Control*. Elsevier, New York.
- Shebuski, J.R., O. Vilhelmsson, and K.J. Miller. 2000. Effects of growth at low water activity on the thermal tolerance of *Staphylococcus aureus*. *Journal of Food Protection* 63:1277-1281.
- Taoukis, P., W. Breene, and T.P. Labuza. (1988). Intermediate moisture foods. *Advances in Cereal Science and Technology*. 9:91-128.
- Tapia de Daza, M.S., Y. Villegas, and A. Martinez. 1991. Minimal water activity for growth of *Listeria monocytogenes* as affected by solute and temperature. *International Journal of Food Microbiology* 14:333-337.
- Tokuoka, K., and T. Ishitani. 1991. Minimum water activities for the growth of yeasts isolated from high-sugar foods. *Journal of General and Applied Microbiology* 37:111-119.
- Torres, R., J. Usall, N. Teixido, M. Abadias, and I. Vinas. 2003. Liquid formulation of the biocontrol agent *Candida sake* by modifying water activity or adding protectants. *Journal of Applied Microbiology* 94:330-339.
- Ucar, F., and I. Guneri. 1996. The effect of water activity, pH and temperature on the growth of osmophilic yeasts. *Turkish Journal of Biology* 20:37-46.
- Wijtzes, T., P.J. McClure, M.H. Zwietering, and T.A. Roberts. 1993. Modelling bacterial growth of *Listeria monocytogenes* as a function of water activity, pH and temperature. *International Journal of Food Microbiology* 18:139-149.
- Zwietering, M.H., T. Wijtzes, J.C. de Wit, and K. Van't Riet. 1992. A decision support system for prediction of the microbial spoilage in foods. *Journal of Food Protection* 55:973-979.

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- Beuchat, L. R. (1981). Microbial stability as affected by water activity. *Cereal Foods World*. 26(7):345- 349.
- Chen, H.C. (1995). Seafood microorganisms and seafood safety. *Journal of Food and Drug Analysis*. 3:133-144.
- Farber, J.M., F. Coates, and E. Daley. (1992). Minimum water activity requirements for the growth of *Listeria monocytogenes*. *Letters In Applied Microbiology*. 15:103-105.
- Garcia de Fernando, G.D., O. Diaz, M. Fernandez, and J.A. Ordonez. (1992). Changes in water activity of selected solid culture media throughout incubation. *Food Microbiology*. 9:77-82.

- Gibson, A.M., J. Baranyi, J.I. Pitt, M.J. Eyles, and T.A. Roberts. (1994). Predicting fungal growth: The effect of water activity on *Aspergillus flavus* and related species. *International Journal of Food Microbiology*. 23:419-431.
- Goalen, N., J.E. Smith, J. Lacey, and G. Gettinby. (1997). Effects of temperature, water activity, and incubation time on production of aflatoxins and cyclopiazonic acid by an isolate of *Aspergillus flavus* in surface agar culture. *Applied and Environmental Microbiology*. 63:1048-1053.
- Hocking, A.D. and B.F. Miscamble. (1995). Water relations of some Zygomycetes isolated from food. *Mycological Research*. 99:1113-1118.
- Hocking, A.D., B.F. Miscamble, and J.I. Pitt. (1994). Water relations of *Alternaria alternata*, *Cladosporium cladosporioides*, *Cladosporium sphaerospermum*, *Curvularia lunata* and *Curvularia pallescens*. *Mycological Research*. 98:91- 94.
- Houtsma, P.C., A. Heuvelink, J. Dufrenne, and S. Notermans. (1994). Effect of sodium lactate on toxin production, spore germination and heat resistance of proteolytic *Clostridium botulinum* strains. *Journal of Food Protection*. 57:327-330.
- Kuntz, L.A. (1992). Keeping microorganisms in control. *Food Product Design*. August:44-51.
- Li, K.Y. and J.A. Torres. (1993). Water activity relationships for selected mesophiles and psychrotrophs at refrigeration temperature. *Journal of Food Protection*. 56:612-615.
- Marauska, M., A. Vigants, A. Klincare, D. Upite, E. Kaminska, and M. Bekers. (1996). Influence of water activity and medium osmolality on the growth and acid production of *Lactobacillus casei* var. *alactosus*. *Proceedings of the Latvian Academy of Sciences Section B Natural Exact and Applied Sciences*. 50:144-146.
- Miller, A.J. (1992). Combined water activity and solute effects on growth and survival of *Listeria monocytogenes* Scott A. *Journal of Food Protection*. 55:414-418.
- Nakajo, M. and Y. Moriyama. (1993). Effect of pH and water activity on heat resistance of spores of *Bacillus coagulans*. *Journal of the Japanese Society for Food Science and Technology*. 40:268-271.
- Nolan, D.A., D.C. Chamblin, and J.A. Troller. (1992). Minimal water activity levels for growth and survival of *Listeria monocytogenes* and *Listeria innocua*. *International Journal of Food Microbiology*. 16:323-335.
- Petersson, S. and J. Schnuerer. (1995). Biocontrol of mold growth in high-moisture wheat stored under airtight conditions by *Pichia anomala*, *Pichia guilliermondii*, and *Saccharomyces cerevisiae*. *Applied and Environmental Microbiology*. 61:1027-1032.
- Pitt, J.I. and B.F. Miscamble. (1995). Water relations of *Aspergillus flavus* and closely related species. *Journal of Food Protection*. 58:86-90.
- Quintavalla, S. and G. Parolari. (1993). Effects of temperature, a w and pH on the growth of *Bacillus* cells and spore: A response surface methodology study. *International Journal of Food Microbiology*. 19:207-216.
- Saad, R.R. (1992). Effect of water activity on growth and lipids of xerophilic fungi, *Aspergillus repens* and *Aspergillus amstelodami*. *Zentralblatt Fuer Mikrobiologie*. 147:61-64.

- Santos, J., T.M. Lopez-Diaz, M.L. Garcia-Lopez, M.C. Garcia-Fernandez, and A. Otero. (1994). Minimum water activity for the growth of *Aeromonas hydrophila* as affected by strain, temperature and humectant. *Letters In Applied Microbiology*. 19:76-78.
- Tapia de Daza, M.S., Y. Villegas, and A. Martinez. (1991). Minimal water activity for growth of *Listeria monocytogenes* as affected by solute and temperature. *International Journal of Food Microbiology*. 14:333-338.
- Tokuoka, K. and T. Ishitani. (1991). Minimum water activities for the growth of yeasts isolated from high sugar foods. *Journal of General and Applied Microbiology*. 37:111-119.
- Ucar, F. and I. Guneri. (1996). The effect of water activity (aw), pH and temperature on the growth of osmophilic yeasts. *Turkish Journal of Biology*. 20:37-46.
- Wijtzes, T., P.J. McClure, M.H. Zwietering, and T.A. Roberts. (1993). Modelling bacterial growth of *Listeria monocytogenes* as a function of water activity, pH and temperature. *International Journal of Food Microbiology*. 18:139-149.
- Zwietering, M.H., T. Wijtzes, J.C. De-Wit, and R.K. Van'T. (1992). A decision support system for prediction of the microbial spoilage in foods. *Journal of Food Protection*. 55:973-979.

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- Allen, K., D. Cornforth, D. Whittier, M. Vasavada, and B. Nummer. 2007. Evaluation of high humidity and wet marinade methods for pasteurization of jerky. *Journal of Food Science*. 72:C351-C355.
- Chen, N. and L.A. Shelef. (1992). Relationship between water activity, salts of lactic acid, and growth of *Listeria monocytogenes* in a meat model system. *Journal of Food Protection*. 55:574-578.
- Clavero, M.R.S. and L.R. Beuchat. (1996). Survival of *Escherichia coli* O157:H7 in broth and processed salami as influenced by pH, water activity, and temperature and suitability of media for its recovery. *Applied and Environmental Microbiology*. 62:2735-2740.
- Duffy, L.L., P.B. Vanderlinde, and F.H. Grau. (1994). Growth of *Listeria monocytogenes* on vacuum-packed cooked meats: Effects of pH, a_w, nitrite and ascorbate. *International Journal of Food Microbiology*. 23:377-390.
- Elgasim, E.A., and M.S. Al Wesali. 2000. Water activity and Hunter colour values of beef patties extended with samh (*Mesembryanthemum forsskalei* Hochst) flour. *Food Chem* 69(2):181-185.
- Fernandez-Salguero J., R. Gómez, and M.A. Carmona. (1994). Water activity of Spanish intermediate- moisture meat products. *Meat Science*. 38:341- 346.
- Gómez, R. and Fernandez-Salguero J. (1993). Note: Water activity of Spanish intermediate moisture fish products. *Revista Espanola De Ciencia Y Tecnologia De Alimentos*. 33:651-656.
- Hand, L. (1994). Controlling water activity and pH in snack sticks. *Meat Marketing and Technology*. May:55-56.

- Lee, M.B. and S. Styliadis. (1996). A survey of pH and water activity levels in processed salamis and sausages in Metro Toronto. *Journal of Food Protection*. 59:1007-1010.
- Luecke, F.K. (1994). Fermented meat products. *Food Research International*. 27:299-307.
- Minegishi, Y., Y. Tsukamasa, K. Miake, T. Shimasaki, C. Imai, M. Sugiyama, and H. Shinano. (1995). Water activity and microflora in commercial vacuum-packed smoked salmons. *Journal of the Food Hygienic Society of Japan*. 36:442-446.
- Nunez, F., M.C. Diaz, M. Rodriguez, E. Aranda, A. Martin, and M.A. Asensio. 2000. Effects of substrate, water activity, and temperature on growth and verrucosidin production by *Penicillium polonicum* isolated from dry-cured ham. *Journal of Food Protection* 63:231-236.
- Placido, M. and M.P. Aleman. 2002. Rapid hygrometric method for determining water activity. *Ciencia y Tecnologia Alimentaria* 3(4):229-235.
- Rocha-Garza, A.E. and J.F. Zayas. (1996). Quality of broiled beef patties supplemented with wheat germ protein flour. *Journal of Food Science*. 61:418-421.
- Sabadini, E., M.D. Hubinger, P.-J.d.Sobral, and B.C. Carvalho, Jr. 2001. Change of water activity and meat colour in the elaboration process of dehydrated salted meat. *Ciencia e Tecnologia de Alimentos* 21(1):14-19.
- Shimasaki, T., K. Miake, Y. Tsukamasa, M.A. Sugiyama, Y. Minegishi, and H. Shinano. (1994). Effect of Water Activity and Storage Temperature on the Quality and Microflora of Smoked Salmon. *Nippon Suisan Gakkaishi*. 60:569-576.
- Untermann, F. and C. Muller. (1992). Influence of a w value and storage temperature on the multiplication and enterotoxin formation of staphylococci in dry-cured raw hams. *International Journal of Food Microbiology*. 16:109-115.
- Williams, S.K., G.E. Rodrick, and R.L. West. (1995). Sodium lactate affects shelf life and consumer acceptance of fresh (*Ictalurus nebulosus*, *marmoratus*) fillets under simulated retail conditions. *Journal of Food Science*. 60:636-639.

13.4.2. Mlékárenské produkty

- Clavero, M.R.S., and L.R. Beuchat. 1996. Survival of *Escherichia coli* O157:H7 in broth and processed salami as influenced by pH, water activity, and temperature and suitability of media for its recovery. *Appl Environ Microbiol* 62:2735-2740.
- Correia, R., M. Magalhaes, M. Pedrini, A. da Cruz, and I. Clementino. 2008. Ice cream made from cow and goat milk: chemical composition and melting point characteristics. *Revista Ciencia Agronomica* 39:251-256.
- Duffy, L.L., P.B. Vanderlinde, and F.H. Grau. 1994. Growth of *Listeria monocytogenes* on vacuum-packed cooked meats: Effects of pH, aw, nitrite and ascorbate. *International Journal of Food Microbiology* 23:377-390.
- Fresno, J.M., M.E. Tornadijo, J. Carballo, P.J. Gonzalez, and A. Bernardo. (1996). Characterization and biochemical changes during the ripening of a Spanish craft goat's milk cheese (Armada variety). *Food Chemistry*. 55:225-230.
- Gomez, R., and J. Fernandez-Salguero. 1993. Note: Water activity of Spanish intermediate moisture fish products. *Revista Espanola De Ciencia Y Tecnologia De Alimentos* 33:651-656.

-
- Hand, L. 1994. Controlling water activity and pH in snack sticks. *Meat Marketing and Technology* May:55-56.
- Hardy, J., J. Scher, and S. Banon. 2002. Water activity and hydration of dairy powders. *Lait* 82:441-442.
- Hong, Y.H. (1991). Physical and chemical properties of the process cheese on the domestic market. *Korean Journal of Animal Science*. 33:387-391.
- Kombila, M.E. and C. Lacroix. (1991). The effect of combinations of salt, lactose and glycerol on the water activity (A W) of cheese spreads. *Canadian Institute of Food Science and Technology Journal*. 24:233-238.
- Lee, M.B., and S. Styliadis. 1996. A survey of pH and water activity levels in processed salamis and sausages in Metro Toronto. *Journal of Food Protection* 59:1007-1010.
- Luecke, F.K. 1994. Fermented meat products. *Food Res Intl* 27:299-307.
- Malec, L.S., A.S. Pereyra-Gonzales, G.B. Naranjo, and M.S. Vigo. 2002. Influence of water activity and storage temperature on lysine availability of a milk like system. *Food Res Intl* 35(9):849-853.
- Minegishi, Y., Y. Tsukamasa, K. Miake, T. Shimasaki, C. Imai, M. Sugiyama, and H. Shinano. 1995. Water activity and microflora in commercial vacuum-packed smoked salmons. *Journal of the Food Hygienic Society of Japan* 36:442-446.
- Pisecky, J. (1992). Water activity of milk powders. *Milchwissenschaft*. 47:3-7.
- Rocha-Garza, A.E., and J.F. Zayas. 1996. Quality of broiled beef patties supplemented with wheat germ protein flour. *Journal of Food Science* 61:418-421.
- Shah, N.P., and R.R. Ravula. 2000. Influence of water activity on fermentation, organic acids production and viability of yoghurt and probiotic bacteria. *Australian Journal of Dairy Technology* 55(3):127-131.
- Shimasaki, T., K. Miake, Y. Tsukamasa, M.A. Sugiyama, Y. Minegishi, and H. Shinano. 1994. Effect of water activity and storage temperature on the quality and microflora of smoked salmon. *Nippon Suisan Gakkaishi* 60:569-576.
- Tornadijo, E., J.M. Fresno, J. Carballo, and S.R. Martin. (1993). Study of Enterobacteriaceae throughout the manufacturing and ripening of hard goats' cheese. *Journal of Applied Bacteriology*. 75:240-246.
- Untermann, F., and C. Muller. 1992. Influence of aw value and storage temperature on the multiplication and enterotoxin formation of staphylococci in dry-cured raw hams. *International Journal of Food Microbiology* 16:109-115.
- Valik, L. and F. Gerner. (1995). Effect of water activity adjusted with different solutes on growth and Lactic acid production by *Lactobacillus helveticus*. *Folia Microbiologica*. 40:472-474.
- Vivier, D., M. Rivemale, J.P. Reverbel, R. Ratomahenina, and P. Galzy. (1994). Study of the growth of yeasts from feta cheese. *International Journal of Food Microbiology*. 22:207-215.
- Vivier, D., R. Ratomahenina, and P. Galzy. (1994). Characteristics of micrococci from the surface of Roquefort cheese. *Journal of Applied Bacteriology*. 76:546-552.
- Williams, S.K., G.E. Rodrick, and R.L. West. 1995. Sodium lactate affects shelf life and consumer acceptance of fresh Catfish (*Ictalurus nebulosus*, marmoratus) fillets under simulated retail conditions. *Journal of Food Science* 60:636-639.
-

13.4.3. Ovoce a zelenina

- Ayub, M., R. Khan, S. Wahab, A. Zeb, and J. Muhammad. (1995). Effect of crystalline sweeteners on the water activity and shelf stability of osmotically dehydrated guava. *Sarhad Journal of Agriculture*. 11:755-761.
- Beveridge, T. and S.E. Weintraub. (1995). Effect of blanching pretreatment on color and texture of apple slices at various water activities. *Food Research International*. 28:83-86.
- Clavero, M.R.S., R.E. Brackett, L.R. Beuchat, and M.P. Doyle. 2000. Influence of water activity and storage conditions on survival and growth of proteolytic *Clostridium botulinum* in peanut spread. *Food Microbiology* 17(1):53-61.
- Fouskaki, M., K. Karametsi, and N.A. Chaniotakis. 2003. Method for the determination of water content in sultana raisins using a water activity probe. *Food Chem* 82:133-1337.
- Gogus, F., C. Cuzdemir, and S. Eren. 2000. Effects of some hydrocolloids and water activity on nonenzymic browning of concentrated orange juice. *Nahrung* 44(6):438-442.
- Hubinger, M., F.C. Menegalli, R.J. Aguerre, and C. Suarez. (1992). Water vapor adsorption isotherms of guava, mango and pineapple. *Journal of Food Science*. 57:1405-1407.
- Jimenez, M., M. Manez, and E. Hernandez. (1996). Influence of water activity and temperature on the production of zearalenone in corn by three *Fusarium* species. *International Journal of Food Microbiology*. 29:417-421.
- Khalloufi, S., J. Giasson, and C. Ratti. 2000. Water activity of freeze dried mushrooms and berries. *Canadian Agricultural Engineering* 42(1):51-56.
- Kiranoudis, C.T., Z.B. Maroulis, E. Tsami, and K.D. Marinos. (1993). Equilibrium moisture content and heat of desorption of some vegetables. *Journal of Food Engineering*. 20:55-74.
- Lopez-Malo, A., and E. Palou. 2000. Modeling the growth/nogrowth interface of *Zygosaccharomyces bailii* in Mango puree. *Journal of Food Science*: 65:516-520.
- Makower, B. and G.L. Dehority. (1943). Equilibrium moisture content of dehydrated vegetables. *Industrial and Engineering Chemistry*. 35(2):193-197.
- Maltini, E., D. Torreggiani, B.R. Brovotto, and G. Bertolo. (1993). Functional properties of reduced moisture fruits as ingredients in food systems. *Food Research International*. 26:413-419.
- Marin, S., N. Magan, M. Abellana, R. Canela, A.J. Ramos, and V. Sanchis. 2000. Selective effect of propionates and water activity on maize mycoflora and impact on fumonisin B1 accumulation. *Journal of Stored Products Research* 36:203-214.
- Marin, S., V. Sanchis, and N. Magan. (1995). Water activity, temperature, and pH effects on growth of *Fusarium moniliforme* and *Fusarium proliferatum* isolates from maize. *Canadian Journal of Microbiology*. 41:1063-1070.
- Monsalve-Gonzalez, A., G.V. Barbosa-Canovas, and R.P. Cavalieri. (1993). Mass transfer and textural changes during processing of apples by combined methods. *Journal of Food Science*. 58:1118-1124.
- Pinsirodom, P., and K.L. Parkin. 2000. Selectivity of Celite-immobilized patatin (lipid acyl hydrolase) from potato (*Solanum tuberosum* L.) tubers in esterification reactions

as influenced by water activity and glycerol analogues as alcohol acceptors. *J. Agric. Food Chem.* 48(2):155-160.

Tapia de Daza, M.S., C.E. Aguilar, V. Roa, and R.V. Diaz de Tablante. (1995). Combined stress effects on growth of *Zygosaccharomyces rouxii* from an intermediate moisture papaya product. *Journal of Food Science.* 60:356-359.

Zeb, A., R. Khan, A. Khan, M. Saeed, and S.A. Manan. (1994). Influence of crystalline sucrose and chemical preservatives on the water activity and shelf stability of intermediate banana chips. *Sarhad Journal of Agriculture.* 10:721-726.

Zhang, X.W., X. Liu, D.X. Gu, W. Zhou, R.L. Wang, and P. Liu. (1996). Desorption isotherms of some vegetables. *Journal of the Science of Food and Agriculture.* 70:303-306.

13.4.4. Pečivo a cereálie

Abellana, M., A.J. Ramos, V. Sanchis, and P.V. Nielsen. 2000. Effect of modified atmosphere packaging and water activity on growth of *Eurotium amstelodami*, *E. chevalieri* and *E. herbariorum* on a sponge cake analogue. *Journal of Applied Microbiology* 88:606-616.

Aramouni, F.M., K.K. Kone, J.A. Craig, and D.-Y.C. Fung. (1994). Growth of *Clostridium sporogenes* PA 3679 in home-style canned quick breads. *Journal of Food Protection.* 57:882-886.

Cahagnier, B., L. Lesage, and M.D. Richard. (1993). Mould growth and conidiation in cereal grains as affected by water activity and temperature. *Letters In Applied Microbiology.* 17:7-13.

Clawson, A.R. and A.J. Taylor. (1993). Chemical changes during cooking of wheat. *Food Chemistry.* 47:337-343.

Fleurat-Lessard, F. 2002. Qualitative reasoning and integrated management of the quality of stored grain: a promising new approach. *Journal of Stored Products Research* 38:191-218.

Gómez, R., Fernandez-Salguero J., M.A. Carmona, and D. Sanchez. (1993). Water activity in foods with intermediate moisture levels: Bakery and confectionery products: Miscellany. *Alimentaria.* 30:55-57.

Guynot, M.E., A.J. Ramos, L. Seto, P. Purroy, V. Sanchis, and S. Marin. 2003. Antifungal activity of volatile compounds generated by essential oils against fungi commonly causing deterioration of bakery products.

Harris, M. and M. Peleg. (1996). Patterns of textural changes in brittle cellular cereal foods caused by moisture sorption. *Cereal Chemistry.* 73:225-231

Hope, R., and N. Magan. 2003. Two-dimensional environmental profiles of growth, deoxynivalenol and nivalenol production by *Fusarium culmorum* on wheat-based substrate. *Lett Appl Microbiol* 37:70-74.

Michniewicz, J., C.G. Biliaderis, and W. Bushuk. (1992). Effect of added pentosans on some properties of wheat bread. *Food Chemistry.* 43:251-257.

Moreno-Contreras, M.D., A.J. Martinez-Yepey, and R.R. Martinez. 2000. Determination of deoxynivalenol (DON) in wheat, barley and corn and its relationship with the levels of total molds, *Fusarium* spp., infestation percentage, and water activity. *Archivos Latinoamericanos de Nutricion.* 50(2):183-186.

- Phoungchandang, S., and J.L. Woods. 2000. Moisture diffusion and desorption isotherms for banana. *Journal of Food Science* 65:651-657.
- Ramanathan, S. and S. Cenkowski. (1995). Sorption isotherms of flour and flow behaviour of dough as influenced by flour compaction. *Canadian Agricultural Engineering*. 37:119-124.
- Roessler, P.F. and M.C. Ballenger. (1996). Contamination of an unpreserved semisoft baked cookie with a xerophilic *Aspergillus* species. *Journal of Food Protection*. 59:1055-1060.
- Schebor, C., and J. Chirife. 2000. A survey of water activity and pH values in fresh pasta packed under modified atmosphere manufactured in Argentina and Uruguay. *Journal of Food Protection* 63:965-969.
- Seiler, D.A.L. (1979). The mould-free shelf life of bakery products. *FMBRA Bulletin*. April(2):71-74.
- Sumner, S.S., J.A. Albrecht, and D.L. Peters. (1993). Occurrence of enterotoxigenic strains of *Staphylococcus aureus* and enterotoxin production in bakery products. *Journal of Food Protection*. 56:722-724.
- Tesch, R., M.D. Normand, and M. Peleg. (1996). Comparison of the acoustic and mechanical signatures of two cellular crunchy cereal foods at various water activity levels. *Journal of the Science of Food and Agriculture*. 70:347-354.
- Weegels, P.L., J.A. Verhoek, A.M.G. de Groot, and R.J. Hamer. (1994). Effects of gluten of heating at different moisture contents: I. Changes in functional properties. *Journal of Cereal Science*. 19:31-38.

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- Cardelli, C., and T.P. Labuza. 2001. Application of Weibull Hazard Analysis to the determination of shelf life of roasted and ground coffee. *Lebensm Wiss Technol* 34:273-278.
- Carson, K.J., J.L. Collins, and M.P. Penfield. (1994). Unrefined, dried apple pomace as a potential food ingredient. *Journal of Food Science*. 59:1213-1215.
- Cavia, M.M., M.A. Fernandez-Muino, J.F. Huidobro, and M.T. Sancho. 2004. Correlation between Moisture and Water Activity of Honeys Harvested in Different Years. *Journal of Food Science* 69:C-368-370.
- Durrani, M.J., R. Khan, M. Saeed, and A. Khan. (1992). Development of concentrated beverages from Anna apples with or without added preservatives by controlling activity of water for shelf stability. *Sarhad Journal of Agriculture*. 8:23-28.
- Ferragut, V., J.A. Salazar, and A. Chiralt. (1993). Stability in the conservation of emulsified sauces low in oil content. *Alimentaria*. 30:67-69.
- Gleiter, R.A., H. Horn, and H.-D. Isengard. 2006. Influence of type and state of crystallization on the water activity of honey. *Food Chem* 96:441-445.
- Hajmeer, M.N., F.M. Aramouni, and E.A.E. Boyle. 2000. Shelf-life of lite syrup after opening and storage at room or refrigerated temperature. *Journal of Food Quality* 23:529-540.
- Ibarz, A., J. Pagan, and R. Miguelsanz. (1992). Rheology of clarified fruit juices: II. Blackcurrant juices. *Journal of Food Engineering*. 15:63-74.

- Khalloufi , S., Y. El-Maslouhi, and C. Ratti. 2000. Mathematical model for prediction of glass transition temperature of fruit powders. *Journal of Food Science* 65:842-848.
- Kusumegi, K., T. Takahashi, and M. Miyagi. (1996). Effects of addition of sodium citrate on the pasteurizing conditions in "Tuyu", Japanese noodle soup. *Journal of the Japanese Society for Food Science and Technology*. 43:740-747.
- Perera, C.O. 2005. Selected quality attributes of dried foods. *Drying Technology* 23:717-730.
- Sa, M.M. and A.M. Sereno. (1993). Effect of temperature on sorption isotherms and heats of sorption of quince jam. *International Journal of Food Science and Technology*. 28:241-248.
- Shafi ur-Rahman, M. 2005. Dried food properties: challenges ahead. *Drying Technology* 23:695-715.

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- Ahlneck, C. and G. Zografi . (1990). The molecular basis of moisture effects on the physical and chemical stability of drugs in the solid state. *International Journal of Pharmaceutics*. 62:87-95.
- Bell, L.N., and K.L. White. 2000. Thermal Stability in Solids as Affected by the Glass Transition. *Journal of Food Science* 65:498-501.
- Cochet, N. and A.L. Demain. (1996). Effect of water activity on production of beta-lactam antibiotics by *Streptomyces clavuligerus* in submerged culture. *Journal of Applied Bacteriology*. 80:333-337.
- Costantino, H.R., R. Langer, and A.M. Klibanov. (1994). Solid-phase aggregation of proteins under pharmaceutically relevant conditions. *Journal of Pharmaceutical Sciences*. 83:1662-1669.
- Enigl, D.C. 2001. Pharmaceutical stability testing using water activity. *European Pharmaceutical Review* 6:46-49.
- Enigl, D.C. and K.M. Sorrels. (1997). Water Activity and Self-Preserving Formulas. In: *Preservative-Free and Self-Preserving Cosmetics and Drugs: Principles and Practice*. Kabara, J.J. and D.S. Orth (ed.) Marcel Dekker, pp. 45-73.
- Friedel, R.R. and A.M. Cundell. (1998). The application of water activity measurement to the microbiological attributes testing of nonsterile over-the-counter drug products. *Pharmacopeial Forum*. 24(2):6087-6090.
- Hageman, M.J. (1988). The role of moisture in protein stability. *Drug Development and Industrial Pharmacy*. 14(14):2047-2070.
- Heidemann, D.R. and P.J. Jarosz. (1991). Formulation studies involving moisture uptake in solid dosage forms. *Pharmaceutical Research*. 8(3):292- 297.
- Kontny, M.J. (1988). Distribution of water in solid pharmaceutical systems. *Drug Development and Industrial Pharmacy*. 14(14):1991-2027.
- Sablani, S.S., K. Al-Belushi, I. Al-Marhubi, and R. Al-Belushi. 2007. Evaluating Stability of Vitamin C in Fortified Formula Using Water Activity and Glass Transition. *International Journal of Food Properties* 10:61-71.
- Zografi , G. (1988). States of water associated with solids. *Drug Development and Industrial Pharmacy*. 14(14):1905-1926.

Zografi, G. and M.J. Kontny. (1986). The interactions of water with cellulose- and starch-derived pharmaceutical excipients. *Pharmaceutical Research*. 3(4):187-193.

13.5. Různé

Bell, L.N. and T.P. Labuza. (1992). Compositional influence on the pH of reduced-moisture solutions. *Journal of Food Science*. 57:732-734.

Bell, L.N. and T.P. Labuza. (1994). Influence of the low-moisture state on pH and its implication for reaction kinetics. *Journal of Food Engineering*. 22:291-312.

Bell, L.N. (1995). Kinetics of non-enzymatic browning in amorphous solid systems: Distinguishing the effects of water activity and the glass transition. *Food Research International*. 28:591-597.

Bhandari, B., and I. Bareyre, 2003. Estimation of crystalline phase present in glucose crystal-solution mixture by water activity measurement. *Lebensm Wiss Technol* 36:729-733(5).

Brake, N.C. and O.R. Fennema. (1993). Edible coatings to inhibit lipid migration in a confectionery product. *Journal of Food Science*. 58:1422-1425.

Dole, M. and L. Faller. (1950). Water sorption by synthetic high polymers. *Journal of the American Chemical Society*. 72:414-419.

Fernandez-Salguero J., R. Gómez, and M.A. Carmona. (1993). Water activity in selected high-moisture foods. *Journal of Food Composition and Analysis*. 6:364-369.

Juhan, K., and G.K. Byung. 2000. Lipase-catalyzed synthesis of lysophosphatidylcholine using organic cosolvent for in situ water activity control. *Journal of American Oil Chemists' Society* 77(7):701-797.

Lima, J.R., S.D.S. Campos, and L.-A.G. Goncalves. 2000. Relationship between water activity and texture of roasted and salted cashew kernel. *Journal of Food Science and Technology* 37(5):512-513.

Lomauro, C.J., A.S. Bakshi, and T.P. Labuza. (1985). Evaluation of food moisture sorption isotherm equations. Part I: Fruit, vegetable and meat products. *Lebensmittel-Wissenschaft und-Technologie*. 18:111-117.

Lomauro, C.J., A.S. Bakshi, and T.P. Labuza. (1985). Evaluation of food moisture sorption isotherm equations. Part II: Milk, coffee, tea, nuts, oilseeds, spices and starchy foods. *Lebensmittel-Wissenschaft und-Technologie*. 18:118-124.

Yasuda, H., H.G. Olf, B. Crist, C.E. Lamaze, and A. Peterlin. (1972). Movement of water in homogeneous water-swollen polymers. In: *Water Structure at the Water Polymer Interface*. Jellinek, H.H.G. (ed.) Plenum Press, New York/ London.

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16.4. Literatura

AOAC, Method 978.18D Preparation of Reference Salt Slushes.1995. Official Methods of Analysis of AOAC International.16th Ed. AOAC International, Arlington VA.

Campbell, G.S. and W.H. Gardner.1971. Psychrometric measurement of soil water potential: temperature and bulk density effects. Soil Sci. Soc. Am. Proc. 35:8-12.

Greenspan, L.1977. Humidity fixed points of binary saturated aqueous solutions. J. Res. National Bureau of Stds. A. Physics and Chem. 51A:89-96.

Lang, A.R.G.1967. Osmotic coefficients and water potentials of sodium chloride solutions from 0 to 40 C. Aust. J. Chem. 20:2017-2023.

Robinson, R.A. and R.H. Stokes.1965. Electrolyte Solutions. Butterworths, London.

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© 2008-2018 Qi Analytical s.r.o. Všechna práva vyhrazena.
Název výrobce: METER Group, Inc.
METER Group, Inc.