AQUA LAB

LITE

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9. Literatura

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- Bousquet-Ricard, M., G. Qualyle, T. Pharm, and J.C. Cheftel. (1980). Comparative study of three methods of determining water activity in intermediate moisture foods. *Lebensmittel- Wissenschaft und-Technologie*. 13:169-173.
- Chirife, J., G. Favetto, C. Ferro-Fontan, and S. Resnik. (1983). The water activity of standard saturated salt solutions in the range of intermediate moisture foods. *Lebensmittel-Wissenschaft und-Technologie*. 16:36-38.
- Duckworth, R. (1975). *Water Relations of Foods*. Academic Press, New York.
- Gomez-Diaz, R. (1992). Water activity in foods: Determination methods. *Alimentaria*. 29:77-82. Gómez, R. and Fernandez-Salguero J. (1992). Water activity and chemical composition of some food emulsions. *Food Chemistry*. 45:91-93.
- Gómez, R. and Fernandez-Salguero J. (1992). Water activity and chemical composition of some food emulsions. *Food Chemistry*. 45:91-93.
- Greenspan, L. (1977). Humidity fixed points of binary saturated aqueous solutions. *Journal of Research of the National Bureau of Standards A.Physics and Chemistry*. 81A:89-96.
- Karmas, E. (1981). Measurement of moisture content. *Cereal Foods World*. 26(7):332-334.
- Kitic, D., D.C. Pereira-Jardim, G. Favetto, S.L. Resnik, and J. Chirife. (1986). Theoretical prediction of the water activity of standard saturated salt solutions at various temperatures. *Journal of Food Science*. 51:1037-1042.
- Labuza, T.P. and R. Contreras-Medellin. (1981). Prediction of moisture protection requirements for foods. *Cereal Foods World*. 26(7):335-343.
- Labuza, T.P., K. Acott, S.R. Tatini, and R.Y. Lee. (1976). Water activity determination: A collaborative study of different methods. *Journal of Food Science*. 41:910-917.
- Prior, B.A. (1979). Measurement of water activity in foods: A review. *Journal of Food Protection*. 42(8):668-674.
- Reid, D.S. (1976). Water activity concepts in intermediate moisture foods. In: *Intermediate Moisture Foods*. Davies, R., G.G. Birch, and K.J. Parker (ed.) Applied Science Publishers, London. pp. 54-65.
- Richard, J. and T.P. Labuza. (1990). Rapid determination of the water activity of some reference solutions, culture media and cheese using a new dew point apparatus. *Sciences des Aliments*. 10:57-64.
- Roa, V. and M.S. Tapia de Daza. (1991). Evaluation of water activity measurements with a dew point electronic humidity meter. *Lebensmittel- Wissenschaft und-Technologie*. 24(3):208-213.
- Roos, K.D. (1975). Estimation of water activity in intermediate moisture foods. *Food Technology*. 29:26-30.
- Scott, V.N. and D.T. Bernard. (1983). Influence of temperature on the measurement of water activity of food and salt systems. *Journal of Food Science*. 48:552-554.
- Snively, M.J., J.C. Price, and H.W. Jun. (1990). A comparison of three equilibrium relative humidity measuring devices. *Drug Development and Industrial Pharmacy*. 16(8):1399-1409.

- Stamp, J.A., S. Linscott, C. Lomauro, and T.P. Labuza. (1984). Measurement of water activity of salt solutions and foods by several electronic methods as compared to direct vapor pressure measurement. *Journal of Food Science*. 49:1139- 1142.
- Stoloff, L. (1978). Calibration of water activity measuring instruments and devices: Collaborative study. *Journal of Association of Official Analytical Chemists*. 61:1166-1178.
- Troller, J.A. (1983). Methods to measure water activity. *Journal of Food Protection*. 46:129-134.
- Troller, J.A. and J.H.B. Christian. (1978). *Water Activity and Food*. Academic Press, New York.
- Troller, J.A. and V.N. Scott. (1992). Measurement of water activity (A_w) and acidity. In: *Compendium of Methods for the Microbiological Examination of Foods*. Vanderzant, C. and D.F. Splittstoesser (ed.) American Public Health Association, Washington, D.C. pp. 135-151.
- van den Berg, C. (1985). Water activity. In: *Concentration and Drying of Foods*. MacCarthy, D. (ed.) Elsevier, London. pp. 11-35.
- van den Berg, C. (1991). Food-water relations: Progress and integration, comments and thoughts. In: *Water Relationships in Foods*. Levine, H. and L. Slade (ed.) Plenum Press, New York.
- van den Berg, C. and Bruin. (1981). Water activity and its estimation in food systems: Theoretical aspects. In: *Water Activity: Influences on Food Quality*. Rockland, L.B. and G.F. Stewart (ed.) Academic Press, New York. pp. 1-61.
- Vega, M.H. and G.V. Barbosa-Canovas. (1994). Prediction of water activity in food systems: A review on theoretical models. *Revista Espanola De Ciencia Y Tecnologia De Alimentos*. 34:368- 388.
- Vega, M.H., B. Romanach, and G.V. Barbosa-Canovas. (1994). Prediction of water activity in food systems: A computer program for predicting water activity in multicomponent foods. *Revista Espanola De Ciencia Y Tecnologia De Alimentos*. 34:427-440.
- Vos, P.T. and T.P. Labuza. (1974). Technique for measurements of water activity in the high a_w range. *Journal of Agricultural and Food Chemistry*. 22:326-327.
- Voysey, P. (1993). An evaluation of the AquaLab CX-2 system for measuring water activity. *Digest, Microbiology Section*. 124

9.2. Kvalita a bezpečnost potravín

- Brandt, L. (1996). Bound for success. Controlling water activity gives technologists the edge in developing safe, shelf-stable foods. *Food Formulating*. September:41-48.
- Chirife, J. and B.M. Del-Pilar. (1994). Water activity, glass transition and microbial stability in concentrated/semimoist food systems. *Journal of Food Science*. 59:921-927.
- Chirife, J. and M.P. Buera. (1995). A critical review of some non-equilibrium situations and glass transitions on water activity values of foods in the microbiological growth range. *Journal of Food Engineering*. 25:531-552.
- Chirife, J. and M.P. Buera. (1996). Water activity, water glass dynamics, and the control of microbiological growth in foods. *Critical Reviews in Food Science and Nutrition*. 36(5):465-513.

- Franks, F. (1982). Water activity as a measure of biological viability and quality control. *Cereal Foods World*. 27(9):403-407.
- Franks, F. (1991). Water activity: a credible measure of food safety and quality? *Trends in Food Science and Technology*. March:68-72.
- Hardman, T.M. (1988). *Water and Food Quality*. Elsevier Press, London.
- Kress-Rogers, E. (1993). Food quality measurement. *Food Industry News*. September:23-26.
- Levine, H. and L. Slade. (1991). *Water Relationships in Foods*. Plenum Press, New York.
- Mannheim, C.H., J.X. Liu, and S.G. Gilbert. (1994). Control of water in foods during storage. *Journal of Food Engineering*. 22:509-532.
- McMeekin, T.A. and T. Ross. (1996). Shelf life prediction: Status and future possibilities. *International Journal of Food Microbiology*. 33:65-83.
- Nelson, K.A. and T.P. Labuza. (1994). Water activity and food polymer science: Implications of state on arrhenius and WLF models in predicting shelf life. *Journal of Food Engineering*. 22:271-289.
- Rockland, L.B. and G.F. Stewart. (1981). *Water Activity: Influences on Food Quality*. Academic Press, New York.
- Rockland, L.B. and S.K. Nishi. (1980). Influence of water activity on food product quality and stability. *Food Technology*. 34:42-59.
- Seow, C.C., T.T. Teng, and C.H. Quah. (1988). *Food Preservation by Moisture Control*. Elsevier, New York.
- Taoukis, P., W. Breene, and T.P. Labuza. (1988). Intermediate moisture foods. *Advances in Cereal Science and Technology*. 9:91-128.

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- Beuchat, L. R. (1981). Microbial stability as affected by water activity. *Cereal Foods World*. 26(7):345-349.
- Chen, H.C. (1995). Seafood microorganisms and seafood safety. *Journal of Food and Drug Analysis*. 3:133-144.
- Farber, J.M., F. Coates, and E. Daley. (1992). Minimum water activity requirements for the growth of *Listeria monocytogenes*. *Letters In Applied Microbiology*. 15:103-105.
- Garcia de Fernando, G.D., O. Diaz, M. Fernandez, and J.A. Ordonez. (1992). Changes in water activity of selected solid culture media throughout incubation. *Food Microbiology*. 9:77-82.
- Gibson, A.M., J. Baranyi, J.I. Pitt, M.J. Eyles, and T.A. Roberts. (1994). Predicting fungal growth: The effect of water activity on *Aspergillus flavus* and related species. *International Journal of Food Microbiology*. 23:419-431.
- Goeleni, N., J.E. Smith, J. Lacey, and G. Gettinby. (1997). Effects of temperature, water activity, and incubation time on production of aflatoxins and cyclopiazonic acid by an isolate of *Aspergillus flavus* in surface agar culture. *Applied and Environmental Microbiology*. 63:1048-1053.
- Hocking, A.D. and B.F. Miscamble. (1995). Water relations of some Zygomycetes isolated from food. *Mycological Research*. 99:1113-1118.

- Hocking, A.D., B.F. Miscamble, and J.I. Pitt. (1994). Water relations of *Alternaria alternata*, *Cladosporium cladosporioides*, *Cladosporium sphaerospermum*, *Curvularia lunata* and *Curvularia pallescens*. *Mycological Research*. 98:91- 94.
- Houtsma, P.C., A. Heuvelink, J. Dufrenne, and S. Notermans. (1994). Effect of sodium lactate on toxin production, spore germination and heat resistance of proteolytic *Clostridium botulinum* strains. *Journal of Food Protection*. 57:327-330.
- Kuntz, L.A. (1992). Keeping microorganisms in control. *Food Product Design*. August:44-51.
- Li, K.Y. and J.A. Torres. (1993). Water activity relationships for selected mesophiles and psychrotrophs at refrigeration temperature. *Journal of Food Protection*. 56:612-615.
- Marauska, M., A. Vigants, A. Klincare, D. Upite, E. Kaminska, and M. Bekers. (1996). Influence of water activity and medium osmolality on the growth and acid production of *Lactobacillus casei* var. *alactosus*. *Proceedings of the Latvian Academy of Sciences Section B Natural Exact and Applied Sciences*. 50:144-146.
- Miller, A.J. (1992). Combined water activity and solute effects on growth and survival of *Listeria monocytogenes* Scott A. *Journal of Food Protection*. 55:414-418.
- Nakajo, M. and Y. Moriyama. (1993). Effect of pH and water activity on heat resistance of spores of *Bacillus coagulans*. *Journal of the Japanese Society for Food Science and Technology*. 40:268-271.
- Nolan, D.A., D.C. Chamblin, and J.A. Troller. (1992). Minimal water activity levels for growth and survival of *Listeria monocytogenes* and *Listeria innocua*. *International Journal of Food Microbiology*. 16:323-335.
- Petersson, S. and J. Schnuerer. (1995). Biocontrol of mold growth in high-moisture wheat stored under airtight conditions by *Pichia anomala*, *Pichia guilliermondii*, and *Saccharomyces cerevisiae*. *Applied and Environmental Microbiology*. 61:1027-1032.
- Pitt, J.I. and B.F. Miscamble. (1995). Water relations of *Aspergillus flavus* and closely related species. *Journal of Food Protection*. 58:86-90.
- Quintavalla, S. and G. Parolari. (1993). Effects of temperature, a w and pH on the growth of *Bacillus* cells and spore: A response surface methodology study. *International Journal of Food Microbiology*. 19:207-216.
- Saad, R.R. (1992). Effect of water activity on growth and lipids of xerophilic fungi, *Aspergillus repens* and *Aspergillus amstelodami*. *Zentralblatt Fuer Mikrobiologie*. 147:61-64.
- Santos, J., T.M. Lopez-Diaz, M.L. Garcia-Lopez, M.C. Garcia-Fernandez, and A. Otero. (1994). Minimum water activity for the growth of *Aeromonas hydrophila* as affected by strain, temperature and humectant. *Letters In Applied Microbiology*. 19:76-78.
- Tapia de Daza, M.S., Y. Villegas, and A. Martinez. (1991). Minimal water activity for growth of *Listeria monocytogenes* as affected by solute and temperature. *International Journal of Food Microbiology*. 14:333-338.
- Tokuoka, K. and T. Ishitani. (1991). Minimum water activities for the growth of yeasts isolated from high-sugar foods. *Journal of General and Applied Microbiology*. 37:111-119.
- Ucar, F. and I. Guneri. (1996). The effect of water activity (aw), pH and temperature on the growth of osmophilic yeasts. *Turkish Journal of Biology*. 20:37-46.

Wijtzes, T., P.J. McClure, M.H. Zwietering, and T.A. Roberts. (1993). Modelling bacterial growth of *Listeria monocytogenes* as a function of water activity, pH and temperature. *International Journal of Food Microbiology*. 18:139-149.

Zwietering, M.H., T. Wijtzes, J.C. De-Wit, and R.K. Van'T. (1992). A decision support system for prediction of the microbial spoilage in foods. *Journal of Food Protection*. 55:973-979.

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Chen, N. and L.A. Shelef. (1992). Relationship between water activity, salts of lactic acid, and growth of *Listeria monocytogenes* in a meat model system. *Journal of Food Protection*. 55:574-578.

Clavero, M.R.S. and L.R. Beuchat. (1996). Survival of *Escherichia coli* O157:H7 in broth and processed salami as influenced by pH, water activity, and temperature and suitability of media for its recovery. *Applied and Environmental Microbiology*. 62:2735-2740.

Duffy, L.L., P.B. Vanderlinde, and F.H. Grau. (1994). Growth of *Listeria monocytogenes* on vacuum-packed cooked meats: Effects of pH, a_w , nitrite and ascorbate. *International Journal of Food Microbiology*. 23:377-390.

Fernandez-Salguero J., R. Gómez, and M.A. Carmona. (1994). Water activity of Spanish intermediate- moisture meat products. *Meat Science*. 38:341- 346.

Gómez, R. and Fernandez-Salguero J. (1993). Note: Water activity of Spanish intermediate moisture fish products. *Revista Espanola De Ciencia Y Tecnologia De Alimentos*. 33:651-656.

Hand, L. (1994). Controlling water activity and pH in snack sticks. *Meat Marketing and Technology*. May:55-56.

Lee, M.B. and S. Styliadis. (1996). A survey of pH and water activity levels in processed salamis and sausages in Metro Toronto. *Journal of Food Protection*. 59:1007-1010.

Luecke, F.K. (1994). Fermented meat products. *Food Research International*. 27:299-307.

Minegishi, Y., Y. Tsukamasa, K. Miake, T. Shimasaki, C. Imai, M. Sugiyama, and H. Shinano. (1995). Water activity and microflora in commercial vacuum-packed smoked salmons. *Journal of the Food Hygienic Society of Japan*. 36:442-446.

Rocha-Garza, A.E. and J.F. Zayas. (1996). Quality of broiled beef patties supplemented with wheat germ protein flour. *Journal of Food Science*. 61:418-421.

Shimasaki, T., K. Miake, Y. Tsukamasa, M.A. Sugiyama, Y. Minegishi, and H. Shinano. (1994). Effect of Water Activity and Storage Temperature on the Quality and Microflora of Smoked Salmon. *Nippon Suisan Gakkaishi*. 60:569-576.

Untermann, F. and C. Muller. (1992). Influence of a w value and storage temperature on the multiplication and enterotoxin formation of staphylococci in dry-cured raw hams. *International Journal of Food Microbiology*. 16:109-115.

Williams, S.K., G.E. Rodrick, and R.L. West. (1995). Sodium lactate affects shelf life and consumer acceptance of fresh (*Ictalurus nebulosus*, *marmoratus*) fillets under simulated retail conditions. *Journal of Food Science*. 60:636-639.

9.4.2. Mlékárenské produkty

- Fresno, J.M., M.E. Tornadijo, J. Carballo, P.J. Gonzalez, and A. Bernardo. (1996). Characterization and biochemical changes during the ripening of a Spanish craft goat's milk cheese (Armada variety). *Food Chemistry*. 55:225-230.
- Hong, Y.H. (1991). Physical and chemical properties of the process cheese on the domestic market. *Korean Journal of Animal Science*. 33:387-391.
- Kombila, M.E. and C. Lacroix. (1991). The effect of combinations of salt, lactose and glycerol on the water activity (A W) of cheese spreads. *Canadian Institute of Food Science and Technology Journal*. 24:233-238.
- Pisecky, J. (1992). Water activity of milk powders. *Milchwissenschaft*. 47:3-7.
- Tornadijo, E., J.M. Fresno, J. Carballo, and S.R. Martin. (1993). Study of Enterobacteriaceae throughout the manufacturing and ripening of hard goats' cheese. *Journal of Applied Bacteriology*. 75:240-246.
- Valik, L. and F. Gerner. (1995). Effect of water activity adjusted with different solutes on growth and Lactic acid production by *Lactobacillus helveticus*. *Folia Microbiologica*. 40:472-474.
- Vivier, D., M. Rivemale, J.P. Reverbel, R. Ratomahenina, and P. Galzy. (1994). Study of the growth of yeasts from feta cheese. *International Journal of Food Microbiology*. 22:207-215.
- Vivier, D., R. Ratomahenina, and P. Galzy. (1994). Characteristics of micrococci from the surface of Roquefort cheese. *Journal of Applied Bacteriology*. 76:546-552.

9.4.3. Ovoce a zelenina

- Ayub, M., R. Khan, S. Wahab, A. Zeb, and J. Muhammad. (1995). Effect of crystalline sweeteners on the water activity and shelf stability of osmotically dehydrated guava. *Sarhad Journal of Agriculture*. 11:755-761.
- Beveridge, T. and S.E. Weintraub. (1995). Effect of blanching pretreatment on color and texture of apple slices at various water activities. *Food Research International*. 28:83-86.
- Hubinger, M., F.C. Menegalli, R.J. Aguerre, and C. Suarez. (1992). Water vapor adsorption isotherms of guava, mango and pineapple. *Journal of Food Science*. 57:1405-1407.
- Jimenez, M., M. Manez, and E. Hernandez. (1996). Influence of water activity and temperature on the production of zearalenone in corn by three *Fusarium* species. *International Journal of Food Microbiology*. 29:417-421.
- Kiranoudis, C.T., Z.B. Maroulis, E. Tsami, and K.D. Marinos. (1993). Equilibrium moisture content and heat of desorption of some vegetables. *Journal of Food Engineering*. 20:55-74.
- Makower, B. and G.L. Dehority. (1943). Equilibrium moisture content of dehydrated vegetables. *Industrial and Engineering Chemistry*. 35(2):193-197.
- Maltini, E., D. Torreggiani, B.R. Brovotto, and G. Bertolo. (1993). Functional properties of reduced moisture fruits as ingredients in food systems. *Food Research International*. 26:413-419.
- Marin, S., V. Sanchis, and N. Magan. (1995). Water activity, temperature, and pH effects on growth of *Fusarium moniliforme* and *Fusarium proliferatum* isolates from maize. *Canadian Journal of Microbiology*. 41:1063-1070.

- Monsalve-Gonzalez, A., G.V. Barbosa-Canovas, and R.P. Cavalieri. (1993). Mass transfer and textural changes during processing of apples by combined methods. *Journal of Food Science*. 58:1118-1124.
- Tapia de Daza, M.S., C.E. Aguilar, V. Roa, and R.V. Diaz de Tablante. (1995). Combined stress effects on growth of *Zygosaccharomyces rouxii* from an intermediate moisture papaya product. *Journal of Food Science*. 60:356-359.
- Zeb, A., R. Khan, A. Khan, M. Saeed, and S.A. Manan. (1994). Influence of crystalline sucrose and chemical preservatives on the water activity and shelf stability of intermediate banana chips. *Sarhad Journal of Agriculture*. 10:721-726.
- Zhang, X.W., X. Liu, D.X. Gu, W. Zhou, R.L. Wang, and P. Liu. (1996). Desorption isotherms of some vegetables. *Journal of the Science of Food and Agriculture*. 70:303-306.

9.4.4. Pečivo a cereálie

- Aramouni, F.M., K.K. Kone, J.A. Craig, and D.-Y.C. Fung. (1994). Growth of *Clostridium sporogenes* PA 3679 in home-style canned quick breads. *Journal of Food Protection*. 57:882-886.
- Cahagnier, B., L. Lesage, and M.D. Richard. (1993). Mould growth and conidiation in cereal grains as affected by water activity and temperature. *Letters In Applied Microbiology*. 17:7-13.
- Clawson, A.R. and A.J. Taylor. (1993). Chemical changes during cooking of wheat. *Food Chemistry*. 47:337-343.
- Gómez, R., Fernandez-Salguero J., M.A. Carmona, and D. Sanchez. (1993). Water activity in foods with intermediate moisture levels: Bakery and confectionery products: *Miscellany. Alimentaria*. 30:55-57.
- Harris, M. and M. Peleg. (1996). Patterns of textural changes in brittle cellular cereal foods caused by moisture sorption. *Cereal Chemistry*. 73:225-231
- Michniewicz, J., C.G. Biliaderis, and W. Bushuk. (1992). Effect of added pentosans on some properties of wheat bread. *Food Chemistry*. 43:251-257.
- Ramanathan, S. and S. Cenkowski. (1995). Sorption isotherms of flour and flow behaviour of dough as influenced by flour compaction. *Canadian Agricultural Engineering*. 37:119-124.
- Roessler, P.F. and M.C. Ballenger. (1996). Contamination of an unpreserved semisoft baked cookie with a xerophilic *Aspergillus* species. *Journal of Food Protection*. 59:1055-1060.
- Seiler, D.A.L. (1979). The mould-free shelf life of bakery products. *FMBRA Bulletin*. April(2):71-74.
- Sumner, S.S., J.A. Albrecht, and D.L. Peters. (1993). Occurrence of enterotoxigenic strains of *Staphylococcus aureus* and enterotoxin production in bakery products. *Journal of Food Protection*. 56:722-724.
- Tesch, R., M.D. Normand, and M. Peleg. (1996). Comparison of the acoustic and mechanical signatures of two cellular crunchy cereal foods at various water activity levels. *Journal of the Science of Food and Agriculture*. 70:347-354.
- Weegels, P.L., J.A. Verhoek, A.M.G. de Groot, and R.J. Hamer. (1994). Effects of gluten of heating at different moisture contents: I. Changes in functional properties. *Journal of Cereal Science*. 19:31-38.

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- Carson, K.J., J.L. Collins, and M.P. Penfield. (1994). Unrefined, dried apple pomace as a potential food ingredient. *Journal of Food Science*. 59:1213-1215.
- Durrani, M.J., R. Khan, M. Saeed, and A. Khan. (1992). Development of concentrated beverages from Anna apples with or without added preservatives by controlling activity of water for shelf stability. *Sarhad Journal of Agriculture*. 8:23-28.
- Ferragut, V., J.A. Salazar, and A. Chiralt. (1993). Stability in the conservation of emulsified sauces low in oil content. *Alimentaria*. 30:67-69.
- Ibarz, A., J. Pagan, and R. Miguelsanz. (1992). Rheology of clarified fruit juices: II. Blackcurrant juices. *Journal of Food Engineering*. 15:63-74.
- Kusumegi, K., T. Takahashi, and M. Miyagi. (1996). Effects of addition of sodium citrate on the pasteurizing conditions in "Tuyu", Japanese noodle soup. *Journal of the Japanese Society for Food Science and Technology*. 43:740-747.
- Sa, M.M. and A.M. Sereno. (1993). Effect of temperature on sorption isotherms and heats of sorption of quince jam. *International Journal of Food Science and Technology*. 28:241-248.

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- Ahneck, C. and G. Zografi . (1990). The molecular basis of moisture effects on the physical and chemical stability of drugs in the solid state. *International Journal of Pharmaceutics*. 62:87-95.
- Cochet, N. and A.L. Demain. (1996). Effect of water activity on production of beta-lactam antibiotics by *Streptomyces clavuligerus* in submerged culture. *Journal of Applied Bacteriology*. 80:333-337.
- Costantino, H.R., R. Langer, and A.M. Klibanov. (1994). Solid-phase aggregation of proteins under pharmaceutically relevant conditions. *Journal of Pharmaceutical Sciences*. 83:1662-1669.
- Enigl, D.C. and K.M. Sorrels. (1997). Water Activity and Self-Preserving Formulas. In: *Preservative-Free and Self-Preserving Cosmetics and Drugs: Principles and Practice*. Kabara, J.J. and D.S. Orth (ed.) Marcel Dekker, pp. 45-73.
- Friedel, R.R. and A.M. Cundell. (1998). The application of water activity measurement to the microbiological attributes testing of nonsterile over-the-counter drug products. *Pharmaceutical Forum*. 24(2):6087-6090.
- Hageman, M.J. (1988). The role of moisture in protein stability. *Drug Development and Industrial Pharmacy*. 14(14):2047-2070.
- Heidemann, D.R. and P.J. Jarosz. (1991). Formulation studies involving moisture uptake in solid dosage forms. *Pharmaceutical Research*. 8(3):292- 297.
- Kontny, M.J. (1988). Distribution of water in solid pharmaceutical systems. *Drug Development and Industrial Pharmacy*. 14(14):1991-2027.
- Zografi , G. (1988). States of water associated with solids. *Drug Development and Industrial Pharmacy*. 14(14):1905-1926.
- Zografi , G. and M.J. Kontny. (1986). The interactions of water with cellulose- and starch-derived pharmaceutical excipients. *Pharmaceutical Research*. 3(4):187-193.

9.5. Různé

- Bell, L.N. and T.P. Labuza. (1992). Compositional influence on the pH of reduced-moisture solutions. *Journal of Food Science*. 57:732-734.
- Bell, L.N. and T.P. Labuza. (1994). Influence of the low- moisture state on pH and its implication for reaction kinetics. *Journal of Food Engineering*. 22:291-312.
- Bell, L.N. (1995). Kinetics of non-enzymatic browning in amorphous solid systems: Distinguishing the effects of water activity and the glass transition. *Food Research International*. 28:591-597.
- Brake, N.C. and O.R. Fennema. (1993). Edible coatings to inhibit lipid migration in a confectionery product. *Journal of Food Science*. 58:1422-1425.
- Dole, M. and L. Faller. (1950). Water sorption by synthetic high polymers. *Journal of the American Chemical Society*. 12:414-419.
- Fernandez-Salguero J., R. Gómez, and M.A. Carmona. (1993). Water activity in selected high-moisture foods. *Journal of Food Composition and Analysis*. 6:364-369.
- Lomauro, C.J., A.S. Bakshi, and T.P. Labuza. (1985). Evaluation of food moisture sorption isotherm equations. Part I: Fruit, vegetable and meat products. *Lebensmittel-Wissenschaft und-Technologie*. 18:111-117.
- Lomauro, C.J., A.S. Bakshi, and T.P. Labuza. (1985). Evaluation of food moisture sorption isotherm equations. Part II: Milk, coffee, tea, nuts, oilseeds, spices and starchy foods. *Lebensmittel- Wissenschaft und-Technologie*. 18:118-124.
- Yasuda, H., H.G. Olf, B. Crist, C.E. Lamaze, and A. Peterlin. (1972). Movement of water in homogeneous water-swollen polymers. In: *Water Structure at the Water Polymer Interface*. Jellinek, H.H.G. (ed.) Plenum Press, New York/ London.

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12.4. Literatura

AOAC, Method 978.18D Preparation of Reference Salt Slushes.1995. Official Methods of Analysis of AOAC International.16th Ed. AOAC International, Arlington VA.

Campbell, G.S. and W.H. Gardner.1971. Psychrometric measurement of soil water potential: temperature and bulk density effects. Soil Sci. Soc. Am. Proc. 35:8-12.

Greenspan, L.1977. Humidity fixed points of binary saturated aqueous solutions. J. Res. National Bureau of Stds. A. Physics and Chem. 51A:89-96.

Lang, A.R.G.1967. Osmotic coefficients and water potentials of sodium chloride solutions from 0 to 40 C. Aust. J. Chem. 20:2017-2023.

Robinson, R.A. and R.H. Stokes.1965. Electrolyte Solutions. Butterworths, London.

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© 2007-2009 Decagon Devices, Inc.
Decagon Devices, Inc